OCULAR INSTRUMENTS

2255 116th Ave NE, Bellevue, Washington 98004-3039 USA

T: 425-455-5200 or 800-888-6616 F: 425-462-6669

E: contact@ocularinc.com I: www.ocularinc.com

Product Care Instructions—English

CLEANING METHOD 1

DEVICE(S): All Ocular Argon/Diode, YAG and Diagnostic Lenses* and OLIV-EQNA, OLIV-WFNA, OLTK-7.2 or 8.2, OTSG, OUV-132-2, OWIV-HMNA, OKSG, OMUSG.

*Note: See Cleaning Method 3 for all O4MAC, Maxfield AC Four Mirror Gonio, and all OG3MAC Autoclavable Three Mirror.

CLEANING:

Rinse: Immediately upon removal from patient's eye, thoroughly rinse in cool or tepid water.

Wash: Place a few drops of mild soap (i.e. neutral ph 7.0) on a moistened cotton ball. Gently clean with a circular motion

Rinse: Thoroughly rinse in cool or tepid water, then dry carefully with a non-linting tissue.

Then: Proceed with either disinfection or sterilization instructions.

DISINFECTION:

Soak in: Glutaraldehyde, 2% or 3.4% aqueous solution, Temperature per manufacturer instructions, Minimum exposure time = 20 minutes OR

Bleach, 10% solution mixed at 1 part bleach to 9 parts water, Recommended exposure time = 10 minutes.

Rinse: Rinse lens thoroughly to remove disinfection solution. 3 cycles of 1 minute, with cool or tepid water is recommended. Dry carefully and place in a dry storage case.

CAUTION: To avoid damage to the lens, do not exceed recommended exposure time.

NOTE: If used on an ulcerated cornea, lens must be STERILIZED before next procedure.

CAUTION: Disinfectant solutions other than those listed above that are approved for use in your area (e.g. FDA, DGHM, VAH, CE Mark...) may be used in accordance with label instructions of the disinfectant manufacturer. Pay strict attention to disinfectant manufacturer recommended concentrations and contact durations. Ensure that disinfectant solution makes complete contact with all device surfaces and lumens.

Note: These lenses are known to be compatible with: Asepti-Wipe, Cavi-cide, Cidex, Cidex OPA, DisCide Wipe, Enviro-cide, Tristel Wipes System, Steris Resert (except OMRA-WFEX, OG4MG-X & OG3MHD-X) and Opti-Cide. Also compatible with H2O2-3%, except the following lenses: OG3M-10, OPDSG, OPDSG-2, OPDSG-3, OS4M, OS4M-2, OK4DG

STERILIZATION:

Autoclave: NO, STERRAD: NO (Caution: Devices will have limited life after 10 cycles), Steris System 1E: Per manufacturer instructions. Steris not compatible with Posners (OPDSG, OPDSG-2, OPDSG-3), Sussmans (OS4M, OS4M-2), OK4DG, OMRA-HM, OMRA-HM-2, OG3M-10, OGP2 & OGP2H. TSO3 Sterizone VP4 125-Z Low Temperature Sterilization System - Cycle 2: Per manufacturer instructions (Caution: Colored aluminum will fade to a natural aluminum color within 25 cycles). EO: Yes - EO Parameters - Minimum Time: 1 hour, Temperature: 130°F (54°C), Aeration Time: 12 hours.

WARNING: Never Steam Autoclave or Boil listed lenses. Never soak in Alcohol, Acetone or Other Solvents.

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