





Ocular Autoclavable Three Mirror Lenses					
	Product Code	Contact OD	Lens Height	Image Mag	Static Gonio FOV
	OG3MAC-10	10mm	25mm	.61x	150°
	OG3MAC-15	15mm	26.5mm		
	OACF-15*	15mm			
	OACF-17*	17mm			
*Flange adapter only   U.S. Patent No. 6,767,098.					
					

### Design

- § The Autoclavable Three Mirror 10mm Lens provides mirrors for the examination of the fundus and the anterior chamber angle.
- § Three mirrors of 64°, 67° and 73° are arranged at 120° intervals.
- § The small 64° mirror is inclined for gonioscopic procedures. It may also be used for the observation of the vitreous and the fundus near the ora serrata.
- § The middle size mirror is inclined at 67° to observe the peripheral fundus from the ora serrata to the region of the equator.
- § The largest mirror is inclined at 73° to observe the fundus from the equator to an area adjacent to the posterior pole.
- § The posterior pole can be observed through the central axis of the lens.
- § A small diameter contact surface for use without methylcellulose.
- § The contact diameter of the OG3MAC-10 can be increased by attaching either the OACF-15(15mm flange diameter) or the OACF-17(17mm flange diameter).

### Cleaning and Disinfection

See Cleaning Method 3

