











Ocular Disposable Vitrectomy Lenses					
Product Code	ODVB	ODVF	ODVM	ODVW	ODV3P
CE 1639 UK CA					
Style	Biconcave	Flat	Magnifying	Wide Field	30° Prism
Image					
Line Drawing					
Static FOV	25°	36°	30°	48°	33°
Image Mag	.80x	1.02x	1.49x	0.49x	1.02x
Description	83D biconcave lens facilitates viewing the fundus in an air-filled vitreous cavity in phakic and pseudophakic eyes.	Plano anterior surface affords a 36 degree field of view of the central posterior pole and vitreous in phakic and pseudophakic eyes. The lens is ideal for photography.	For detailed examination and minute surgical manipulation of retinal membranes in phakic and pseudophakic eyes.	Concave anterior surface facilitates a 48 degree field of view when visualizing the central posterior pole and central vitreous in phakic and pseudophakic eyes.	Provides visualization of the posterior peripheral fundus and vitreous beyond the equator with minimal distortion in phakic, aphakic, and pseudophakic eyes.

Caution: Federal law restricts this device to sale by or on the order of a physician.

Intended Use

Ocular Disposable Vitrectomy lenses are intended to be placed on the eye of a patient to allow a physician to view the posterior structures of the eye during examination and therapy of intraocular abnormalities.

Indication for Use

Ocular Disposable Vitrectomy lenses are indicated for use during vitrectomy and other posterior segment surgeries. The Disposable Vitrectomy lenses provide a clear view of the ocular fundus and retina while viewed through a surgical microscope. Lenses are to be used once and discarded.

Design Features

- Designed to be used once, then discarded.
- Packaged individually in sterile peel pack.
- Silicone flange replaces the need for a suture-down ring. However, the silicone flange can be removed if use of a suture-down ring (such as the Ocular series of lens rings) is preferred.

Method of Use

Rinse the lens with sterile water before use. Place a few drops of an appropriate viscoelastic coupling fluid on the contact surface of the lens prior to placing the lens on the patients eye

Warnings

Do not use if there is damage to lens or sterile package.

Contraindications

The Ocular Disposable Vitrectomy lenses are is intended for transient use only

Cleaning & Sterilization

Lens is provided sterile.

