



Ocular Latina SLT Gonio Laser Lens
Ocular Latina 5 Bar SLT Laser Lens
Ocular Latina 5 Bar Indexing SLT Laser Lens

	Product Code	Style	Gonio Mag	Gonio Laser Spot Mag	Static FOV	Contact OD	Lens Height		
	OLSLT		Flange	1.0x	1.0x	130°	14.5mm		24mm
	OLSLTF						18.0mm		25mm
	OL5SLT		14.5mm				24mm		
	OL5SLTF	Flange	18.0mm				25mm		
	OL5SLT-IR		14.5mm				33mm		
	OL5SLTF-IR	Flange	18.0mm				34mm		

Design

- § The Latina SLT Gonio Laser Lens was designed specifically for Selective Laser Trabeculoplasty.
- § 1.0x magnification maintains laser spot size and 1 to 1 laser energy delivery.
- § Tilted anterior lens surface corrects astigmatism to maintain circular laser beam profile and give sharp images for examination.
- § Suitable for standard laser trabeculoplasty.
- § Large 63° mirror yields bright image for angle photography.
- § OL5SLT & OL5SLTF have five high contrast bars that are imaged over the trabecular meshwork. Bars provide a reference to placement of sub-threshold laser spot.
- § A special lid flange on the OLSLTF & OL5SLTF renders the lens resistant to rejection by the squeezing patient.
- § OL5SLT-IR & OL5SLTF-IR incorporates a hand held Indexing ring that rotates the lens & 5 bars to 10 incremental positions. Provides rotational positioning of 10 laser spots in each of the 10 locations.

Caution

When using the lens for laser photocoagulation, use extreme care to keep the laser away from the edges. If the beam strikes the area around the mirror, it may be absorbed and burn the area. Mirrors damaged in this manner cannot be repaired.

Cleaning & Disinfection

See Cleaning Method 1 included with lens.

- For the OL5SLT-IR & OL5SLTF-IR, remove engraved ring for proper cleaning and disinfection. To remove ring, squeeze two plastic tabs together and slide off. Assemble ring onto lens as pictured.

