
Ocular Brasse Wide Field Capsulotomy Lens

	Product Code	Image Mag	Laser Spot Mag	Contact OD	Lens Height	<i>Designed with:</i> Karl Brasse, MRCOphth, M.D. Muenster, Germany Vienna, Austria Boston, USA	
	OBWAC	1.4x	.70x	16mm	17.2mm		
	UKCA 						

Intended Use

- The Brasse Capsulotomy lens is an ophthalmic contact lens designed to enhance visualization from the anterior chamber to the anterior vitreous.

Indication for Use

- This lens is intended for use by a licensed physician in a manner consistent with standard direct-image ophthalmic contact lens practices.

Lens Design

- The design of the Brasse Capsulotomy lens supports effective reprocessing through steam sterilization.
- A turquoise-colored cone serves as a clear visual indicator to distinguish steam-compatible models from conventional non-autoclavable designs.
- This lens provides an expanded field of view, facilitating comprehensive YAG laser treatment across a larger section of the capsular membrane.

Warnings

- Do not use if the lens is fractured, chipped, scratched, or otherwise damaged.
- The lens must be thoroughly cleaned and sterilized prior to use.

Contra-Indications

- These lenses are designed for temporary application and should not be left on the eye for an extended duration.
- Any lens that touches an ulcerated cornea requires full sterilization before its next application.
- Ensure the correct type and a sufficient volume of coupling fluid is applied between the lens and the cornea before use.

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