
Ocular Ritch Trabeculoplasty Laser Lens

	Product Code	Gonio Mag	Gonio Laser Spot Mag	Contact OD	Lens Height	Static Gonio FOV	<i>Designed with: Robert Ritch, M.D., New York, NY</i> Reference: Review of Ophthalmology Vol. 4, No. 6, p. 97, June 1997	
	ORTA CE	1.4x	.71x	18mm	23mm	80°		


Design

- § The Ritch Trabeculoplasty Laser Lens was designed primarily for anterior chamber photocoagulation and observation procedures.
- § The lens provides four mirrors; two at 59° and two at 64°, plus a central axis view.
- § Positioned over each one of the 59° and 64° mirrors is an anti-reflective coated convex button lens which affords high 1.4x magnification and laser spot size control.
- § The focal distance of the convex button lenses has been calculated to focus strategically on the region of the trabecular meshwork.
- § The lineal height of the lens itself has been shortened so the mirrors are positioned tangent to the anterior lens surface.
- § Adjacent to the mirrors covered by the convex button lenses are two plain mirrors covered with anti-reflective coated glass.
- § The 59° mirror is the optimum mirror for examining and treating the inferior 180° of the angle, while the 64° mirror is best for examining and treating the superior 180° of the angle.
- § While performing laser treatment through the buttons, the spot size is reduced by the effect of the button from 50 microns to 35 microns. The 15 micron spot size reduction increases the laser energy by a factor of two.
- § A combined broad band, anti-reflective coating for both argon and diode lasers is standard.

NOTE

When using lens for photocoagulation, use extreme care to keep the laser beam away from mirror edges. If the beam strikes the black area around the mirror, it can be absorbed and burn the area. Mirrors damaged in this way cannot be repaired.

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© 2001 Ocular Instruments