



## Khaw Surgical Gonioprism

	Product Code	Contact OD	<i>Designed with: P.T Khaw, M.D. PhD London, U.K.</i>	
	OKSG CE	11.5mm		

### Lens Design

- The Khaw surgical gonioprism was designed for easy manipulation during goniotomy and direct viewing gonioscopy procedures.
- A fixation ring attached to the posterior portion of the lens helps to stabilize the globe during surgical procedures.
- Wide field of view lens provides a clear view of the anterior chamber and anterior chamber angle during goniotomy surgical procedures.
- Black anodized aluminum handle provides easy manipulation of the lens.
- Lens magnification provides detailed view of the anterior chamber.
- Use of a coupling viscoelastic may be of considerable help visualizing the anterior chamber angle and to prevent air bubbles in the lens corneal interface.

### Cleaning and Disinfection

See Cleaning Method 8



2255 116th Ave NE, Bellevue, Washington 98004-3039 USA  
T: 425-455-5200 or 800-888-6616 F: 425-462-6669  
E: [contact@ocularinc.com](mailto:contact@ocularinc.com) I: [www.ocularinc.com](http://www.ocularinc.com)

© 2004 OCULAR INSTRUMENTS  
Revised on 2021-06-07 8487C11435