



Ocular Kaufman Two Mirror Laser Lens

	Product Code	Contact OD	Lens Height	
	OK2MA CE	13mm	19.5mm	

Lens Design

- § Two Mirror, with mirrors set at 67° and 76°, is designed for visualization and laser procedures of the retina in all species of monkey.
- § Ocular's Laserlight® high efficiency, broad band, anti-reflective coating provides optimal image contrast, minimizes bothersome reflections and maximizes visible and diode laser transmission.

Cleaning and Disinfection

See Cleaning Method 1



2255 116th Ave NE, Bellevue, Washington 98004-3039 USA
T: 425-455-5200 or 800-888-6616 F: 425-462-6669
E: ocular@ocularinc.com I: www.ocularinc.com

© 2002 Ocular Instruments
7721E3275