



## Ocular Indirect Vitrectomy Lenses

	Product Code	Image Mag	Static FOV	Dynamic FOV	
	OIV-132	.45x	99°	135°	

### Design

Designed to be used on the Ocular Landers Wide Angle Surgical Viewing System (OSVS) in conjunction with an Inverter Vitrectomy System.

Attaches to the OSVS with the Ocular 132D Indirect Vitrectomy Lens Holder (OIV-H132)

Non-contact design allows patient's eye to be rotated freely to view the peripheral retina and vitreous

Steam sterilizable; sterilizing case included

Designed to allow a clear view in the fluid or air filled eye

### Cleaning and Disinfection

See Cleaning Method 3



2255 116th Ave NE, Bellevue, Washington 98004-3039 USA  
 T: 425-455-5200 or 800-888-6616 F: 425-462-6669  
 E: [contact@ocularinc.com](mailto:contact@ocularinc.com) I: [www.ocularinc.com](http://www.ocularinc.com)

© 2001 Ocular Instruments  
 Revised on 2017-08-22 7447H4158