



Ocular 132D Vitrectomy Lens Holders

	Product Code	
	← OIV-H132 OUV-H132-2 →	
	CE	

Design – OIV-H132

Clip style holder for the 132D Indirect Vitrectomy Lens (OIV 132)

Includes two adjustable links that snap onto the end of the slotted arm of the Surgical Viewing System (OSVS)

Steam sterilizable

Design – OUV-H132-2

Holder for the Peyman-Wessels-Landers 132D Upright Vitrectomy Lens (OUV 132-2)

Includes two adjustable links that snap onto the end of the slotted arm of the Surgical Viewing System (OSVS)

Steam sterilizable

Cleaning and Disinfection

See Cleaning Method 3



2255 116th Ave NE, Bellevue, Washington 98004-3039 USA
T: 425-455-5200 or 800-888-6616 F: 425-462-6669
E: contact@ocularinc.com I: www.ocularinc.com

© 2001 Ocular Instruments
Revised on 2017-08-22 7446E4158