Ocular Landers Wide Angle Surgical Viewing System			
	Product Code	Description	Designed With:
	OSVS	Surgical Viewing System with Case	Maurice B. Landers, III, M.D. Houston, TX
Also Available in Sets:			
	OSVS-U132-2	Surgical Viewing System with Ocular Peyman-Wessels-Landers 132D Upright Vitrectomy Lens	
	OSVS-I132	Surgical Viewing System with Ocular 132D Indirect Vitrectomy Lens	C€

Design

Non-contact vitrectomy system designed with a flexible arm for positioning wide angle lenses

Easily swings in and out of surgical field

Can hold indirect lens for use with separate inverter

Designed to be used with the Peyman-Wessels-Landers Wide Angle Lens (OUV 132-2)

More affordable than similar systems

Cleaning and Disinfection

See Cleaning Method 3

