


Ocular Hoskins-Barkan Goniotomy Lenses

	Product Code	Style	Size	<i>Designed with: H. Dunbar Hoskins, M.D San Francisco, CA</i>
	OHBG-1	Infant	10mm	
	OHBG-2	Premature Infant	9mm	
	OHBG-3	Adult	11.5 mm	
CE				

Lens Design

- The Hoskins-Barkan Goniotomy Lens was primarily designed for transverse goniotomy surgery with the operating microscope, but can also be used as a diagnostic lens.
- The infant lens is oval and conical in shape, with a 10mm diameter magnified view of the anterior chamber and anterior chamber angle.
- The premature infant lens is the same in shape and design except the dimensions are 1mm smaller for premature infant surgery.
- There is also an adult size at 11.5 mm diameter.

Cleaning & Disinfection

See Cleaning Method 8



2255 116th Ave NE, Bellevue, Washington 98004-3039 USA
 T: 425-455-5200 or 800-888-6616 F: 425-462-6669
 E: contact@ocularinc.com I: www.ocularinc.com

© 2001 Ocular Instruments
 Revised on 2021-06-07 5505G11432