Ocular Ritch Panoramic Surgical Gonioprism							
T GOND THE	Product Code	Image Mag	Contact OD	Handle Length	Static Field of View		
The state of the s	ORPSG	.73x	10.8mm	77.5mm	160°	Designed with: Robert Ritch, M.D., New York, NY	

Lens Design

- § The Ritch Panoramic Surgical Gonioprism is a glass lens designed for easy manipulation for goniotomy and direct viewing gonioscopy for adults and children.
- § 180° goniotomy can be done.
- § It provides a clear view of the angle during goniotomy surgical procedures.
- § Lens shape creates exposed portion of cornea allowing use of instruments, incisions and corneal retraction sutures.

Cleaning and Disinfection

See Cleaning Method 3

