


## Ocular Reichel OCT SteadyEye

	Product Code	Contact Diameter	Handle Length	<i>Designed with: Elias Reichel, M.D., Boston, MA</i>
	OROCTS	17.5mm	95mm	

### Design

- § Manually stabilizes eye for Optical Coherence Tomography, OCT, in the OR, providing for enhanced image quality.
- § Eliminates eye movement due to microsaccades and/or breathing
- § Sized so it can be used during vitrectomy procedures without interfering with infusion ports
- § Can be used to move the eye in order to image different locations
- § Centers optical axis without interfering with OCT view
- § Functional for children as well as adults

### Cleaning & Disinfection

See Cleaning Method 3



2255 116th Ave NE, Bellevue, Washington 98004-3039 USA  
T: 425-455-5200 or 800-888-6616 F: 425-462-6669  
E: [ocular@ocularinc.com](mailto:ocular@ocularinc.com) I: [www.ocularinc.com](http://www.ocularinc.com)

© 2001 Ocular Instruments  
11355-3671